



INTERNATIONAL CANINE SEMEN BANK - OREGON  
MAIN OFFICE | P.O. Box 651|Sandy, Oregon 97055  
Phone: 503-663-7031 | FAX: (503) 676-8025  
Email: ik9sb@aol.com | Website: www.ik9sb.com

## **ICSB Alkaline pH Buffer**

### **For Adjusting Acidity in Canine Semen Collections**

In cases where acidity is noted of the canine ejaculate, the ICSB Alkaline pH Buffer can protect the sperm cells from the effects of acidity. Semen can be adversely affected in acidic ejaculates, resulting in decreased motility and speed, even to the point of no movement. ICSB Alkaline pH-Adjusting Buffer helps to buffer the acidic pH level less than 6.0 to the optimal pH range of 6.4-6.8 so that the sperm can reach and fertilize the eggs!

ICSB recommends this buffer for use in all acidic ejaculates, whether it will be fresh-collected or fresh-chilled. Keep the buffer frozen until just prior to use. When needed, thaw the buffer at room to body temperature (70-100°F). It is best to add small increments of the buffer, testing the pH throughout with the included pH strips, until the optimal pH level is achieved.

Please contact ICSB if you have any questions or concerns!



INTERNATIONAL CANINE SEMEN BANK - OREGON  
MAIN OFFICE | P.O. Box 651|Sandy, Oregon 97055  
Phone: 503-663-7031 | FAX: (503) 676-8025  
Email: ik9sb@aol.com | Website: www.ik9sb.com

## **ICSB Alkaline pH Buffer**

### **For Adjusting Acidity in Canine Semen Collections**

In cases where acidity is noted of the canine ejaculate, the ICSB Alkaline pH Buffer can protect the sperm cells from the effects of acidity. Semen can be adversely affected in acidic ejaculates, resulting in decreased motility and speed, even to the point of no movement. ICSB Alkaline pH-Adjusting Buffer helps to buffer the acidic pH level less than 6.0 to the optimal pH range of 6.4-6.8 so that the sperm can reach and fertilize the eggs!

ICSB recommends this buffer for use in all acidic ejaculates, whether it will be fresh-collected or fresh-chilled. Keep the buffer frozen until just prior to use. When needed, thaw the buffer at room to body temperature (70-100°F). It is best to add small increments of the buffer, testing the pH throughout with the included pH strips, until the optimal pH level is achieved.

Please contact ICSB if you have any questions or concerns!